

## *InVitroGRO*<sup>™</sup> CP Medium

Product No.	Description	Size
Z99029	<i>InVitroGRO</i> CP Medium	250 mL
Z990003	<i>InVitroGRO</i> CP Medium	500 mL
Z990004	<i>InVitroGRO</i> CP Medium	1 L

### Product Description:

*InVitroGRO* CP Medium is a medium used for plating BioreclamationIVT cryoplateable hepatocytes<sup>1</sup>. *InVitroGRO* CP Medium is formulated without antibiotics to increase its shelf life. BioreclamationIVT recommends the use of the *Torpedo*<sup>™</sup> Antibiotic Mix with *InVitroGRO* CP Medium.

**Stability:** 6 months for *InVitroGRO* CP Medium without *Torpedo* Antibiotic Mix  
7 days for complete *InVitroGRO* CP Medium

**Storage:** 4°C

### Related Products:

Product No.	Size
Z99000	<i>Torpedo</i> Antibiotic Mix 5.5 ml
Z990007	<i>Torpedo</i> Antibiotic Mix 11 mL
Z990008	<i>Torpedo</i> Antibiotic Mix 22 mL

### Reference:

1. Roymans, D.; Van Looveren, C.; Leone, A.; Parker, J. B.; McMillan, M.; Johnson, M. D.; Koganti, A.; Gilissen, R.; Silber, P.; Mannens, G.; Meuldermans, W. Determination of cytochrome P450 1A2 and cytochrome P450 3A4 induction in cryopreserved human hepatocytes. *Biochem. Pharmacol.* **2004**, *67*(3), 427-437.

**Caution:** Treat all products containing human and monkey-derived materials a potentially infectious, as no known test methods can offer assurance that products derived from human or monkey tissues will not transmit infectious agents.

All products are for research use only. Do not use in animals or humans. These products have not been approved for any diagnostic or clinical procedures.