

## *InVitroGRO*<sup>™</sup> Hepatocyte Media

*InVitroGRO*<sup>™</sup> hepatocyte media are used for incubating, thawing, and plating hepatocytes in culture. *InVitroGRO* hepatocyte media are antibiotic-free, giving the media a longer shelf life. *Torpedo*<sup>™</sup> Antibiotic Mix is offered separately, so you can add antibiotics to the *InVitroGRO* hepatocyte media at your convenience. BioreclamationIVT specifically designed the *Torpedo* Antibiotic Mix to complement the *InVitroGRO* hepatocyte media.

### Features and Benefits

- **Optimized.** Media are formulated for culturing hepatocytes
- **Convenient Packaging.** Media have a 6-month shelf life
- **Stability.** Longer shelf life allows you to order in bulk, which saves on shipping costs
- **Ready to Use.** Inventory is immediately available for fast delivery

# InVitroGRO™ Hepatocyte Media

## Ordering Information

Product No.	Description	Use
Z99009	<i>InVitroGRO</i> HI Medium, serum-free, 250 mL	Incubation
Z990012	<i>InVitroGRO</i> HI Medium, serum-free, 500 mL	Incubation
Z990002	<i>InVitroGRO</i> HI Medium, serum-free, 1 L	Incubation
Z99029	<i>InVitroGRO</i> CP Medium, 250 mL	Plating
S01494	<i>InVitroGRO</i> CP Rat Medium, 250 mL	Plating
Z990003	<i>InVitroGRO</i> CP Medium, 500 mL	Plating
Z990004	<i>InVitroGRO</i> CP Medium, 1 L	Plating
Z99019	<i>InVitroGRO</i> HT Medium, 250 mL	Thawing
Z990005	<i>InVitroGRO</i> HT Medium, 500 mL	Thawing
Z990006	<i>InVitroGRO</i> HT Medium, 1 L	Thawing
Z99000	<i>Torpedo</i> Antibiotic Mix 5.5 mL;	Add to <i>InVitroGRO</i> HI, CP, and add to 250 mL of medium HT media
S01495	<i>Torpedo</i> Antibiotic Mix, rat, 5.5 mL; add to 250 mL medium	Add to <i>InVitroGRO</i> CP Rat Medium
Z990007	<i>Torpedo</i> Antibiotic Mix 11 mL;	Add to <i>InVitroGRO</i> HI, CP, and add to 500 mL of medium HT media
Z990008	<i>Torpedo</i> Antibiotic Mix 22 mL;	Add to <i>InVitroGRO</i> HI, CP, and add to 1 L of medium HT media
Z990017	<i>InVitroGRO</i> Hepatocyte Media Kit	Contains one 250-mL bottle each of <i>InVitroGRO</i> HI, CP, and HT media and three 5.5-mL vials of <i>Torpedo</i> Antibiotic Mix
Z99074	<i>InVitroGRO</i> KHB, 250 mL	Buffer

## Technical Information

*InVitroGRO* hepatocyte media are available in 250 mL, 500 mL, and 1 L sizes. *Torpedo* Antibiotic Mix is available for each size of *InVitroGRO* hepatocyte media. *InVitroGRO* hepatocyte media are sterile-filtered and must be stored at 4°C. Once combined with the *Torpedo* Antibiotic Mix, the *InVitroGRO* hepatocyte media have a 7-day shelf life.



*InVitroGRO*™ Hepatocyte Media

*InVitroGRO*™ and *Torpedo*™ are trademarks of BioreclamationIVT Inc.

For Inquiries: North America  
T +1 410 394 7600 F +1 410 455 1245 E [customerservice@bioreclamationivt.com](mailto:customerservice@bioreclamationivt.com)  
Europe  
T +32 2 808 8278 F +32 2 706 5017 E [customerservice@bioreclamationivt.com](mailto:customerservice@bioreclamationivt.com)

 **BioreclamationIVT**  
The complete resource for all biologicals