

InVitroGRO[™] HT Medium

Product No.	Description	Size
Z99019	<i>InVitroGRO</i> HT Medium	250 mL
Z990005	<i>InVitroGRO</i> HT Medium	500 mL
Z990006	<i>InVitroGRO</i> HT Medium	1 L

Product Description:

InVitroGRO HT Medium is a medium used when thawing cryopreserved hepatocytes¹. *InVitroGRO* HT Medium is formulated without antibiotics to increase its shelf life. BioreclamationIVT recommends the use of the *Torpedo*[™] Antibiotic Mix with *InVitroGRO* HT Medium.

Stability: 6 months for *InVitroGRO* HT Medium without *Torpedo* Antibiotic Mix
7 days for complete *InVitroGRO* HT Medium

Storage: 4°C

Related Products:

Product No.	Size
Z99000	<i>Torpedo</i> Antibiotic Mix 5.5 ml
Z990007	<i>Torpedo</i> Antibiotic Mix 11 mL
Z990008	<i>Torpedo</i> Antibiotic Mix 22 mL

Reference:

1. Prabhu, S.; Fackett, A.; Lloyd, S.; McClellan, H. A.; Terrell, C. M.; Silber, P. M.; Li, A. P. Identification of glutathione conjugates of troglitazone in human hepatocytes. *Chem. Biol. Interact.* **2002**, 142(1–2), 83–97.

Caution: Treat all products containing human and monkey-derived materials a potentially infectious, as no known test methods can offer assurance that products derived from human or monkey tissues will not transmit infectious agents.

All products are for research use only. Do not use in animals or humans. These products have not been approved for any diagnostic or clinical procedures.